## Laser Beam Stabilization System

Laser-Laboratorium Göttingen e.V.

The limited pointing stability of high power laser systems is a question of major concern, since fluctuations of the lateral or angular beam position can cause tremendous problems especially in industrial applications. For compensation of drift effects a Laser Beam Stabilization System was developed, which may be utilized for stabilization and realignment of practically all lasers operating in the UV-VIS-NIR spectral range.

## **Characteristics:**

- Piezo driven adaptive mirror
- Closed-loop tip/tilt correction
- Pointing stability: <1 µrad</p>

## **Specifications**

- Accuracy: < 50 mrad</p>
- Max. beam diameter: 40 mm
- Max. divergence : 10 mrad
- Wavelengths from 190 to 1100 nm



The system makes use of a USB camera which allows monitoring of the far-field of the laser beam to be stabilized at a framerate of 400 Hz. For the measurement a small part of the primary beam is coupled into the detection system. The far-field centroid data provide information on the actual beam direction, allowing compensation in x and ydirection by the help of a piezo driven mirror. Automated readjustment after a lateral or angular movement of the beam position is achieved by use of a comprehensive computer program, which allows also determination of the beam pointing stability in accordance with the respective ISO standard 11670.



Laser-Laboratorium Göttingen e.V. Laser-Laboratorium Göttingen e.V. Hans-Adolf-Krebs-Weg 1 D-37077 Göttingen http://www.llg-ev.de/ Dr. Klaus Mann kmann@llg-ev.de Tel: +49-(0)551-5035-41 Fax: +49-(0)551-5035-99